

### KRUUSE Rehab Hock Protector

Designed by ReDog of Sweden™

The hock protector can be used before and after surgery for chronic arthritis i.e. Osteoarthritis, OCD, ligament injuries, arthritis, ED and other assessments.

#### User Manual


The hock protector is placed with the tip of the hock out through the opening in the protector and then stretched around the leg above and under the joint.

#### Care

Wash by hand and air-dry.

#### Size Chart

Size	Cat. No	Suitable for
XXX-Small	279933	Bichon Frisé, Chihuahua, Pomeranian, Maltese
XX-Small	279932	Boston Terrier, Cairn Terrier, Lhasa Apso, Shih Tzu
X-Small	279839	Cocker Spaniel, Whippet, Irish Terrier, Cavalier King Charles Spaniel
Small	279840	Collie, Springer Spaniel, Border Collie, Beagle
Medium	279841	Labrador Retriever, Dalmatian, Golden Retriever, Boxer, Australian Shepherd
Large	279842	Rottweiler, German Shepherd, Golden Retriever, Rhodesian Ridgeback
X-Large	279884	Large Rottweiler, Large German Shepherd, Large Golden Retriever, Great Dane



	Size chart KRUUSE Rehab Hock Protectors in cm	XXXS (cm)	XXS (cm)	XS (cm)	S (cm)	M (cm)	L (cm)	XL (cm)
A	Total height, from edge to edge	8	12	14	16	17.5	19	20.5
B	Total upper width 1 cm from edge	6	8	10	12	14	15	16
C	Total lower width 1 cm from edge	4	5.5	7.5	9.5	11.5	12.5	13.5